



# **Easy 3-Step Operation**


---

**Vet Chemistry Analyzer**

---

## The Easy 3-Step Operation

Dilution, centrifugation, quality control and test result printing will automatically occur after the sample analysis process.



### 1. Insert sample (100 $\mu$ l)

A whole blood sample does not require pre-treatment. You will only need the 100 $\mu$ l of whole blood, serum or plasma to complete the test.



### 2. Insert Rotor (12 minutes running)

This is easy to use in a mobile setting or in the field. The key moving parts are original and sourced from US manufacturers. All fluids are self-contained within the rotor, eliminating the need for internal pumps, valves and other vulnerable parts.



### 3. Read result (Immediately)

The photoelectric colorimetric principle utilized, eliminates random error and cross contamination. The analyzer does quality checks and the calibration occurs automatically and in real time during analysis.



# Contact US



[info@peaclaboratory.co.za](mailto:info@peaclaboratory.co.za)



+27 78 177 7824



[www.peacinfinitum.co.za](http://www.peacinfinitum.co.za)